

VENTON IMPLANT MOTORS

M1

Powerful - Rational - Precise



M3

High-end performance

VENTON M1 LED IMPLANTOLOGY MOTOR

Powerful - Rational - Precise



Warranty :
2 years for the device
1 year for the contra-angle

Maximum torque	70 Ncm
Maximum output power	70 W
Adjustable motor range (rpm)	300 - 40.000 (rpm)
Coolant volume at 100%	> 90ml/min
Transmission ratio	1:1 20:1
Program settings	5 programs for individual adjustment
Main unit dimension	270x240x520mm
Net weight	5.7kg
PV	220V+/-10% - 50~60Hz



Motor

Ease of use, optimization and security
Customizable and free programming to meet different user needs

Foot control

Functions: Water Volume switch, program conversion, motor rotating direction control, stepless speed control




Micro-motor

Compact design, ideal balance and comfort
5-70 Ncm Torque range 5-70 Ncm with 20:1 handpiece

Green ring contra-angle 20:1 LED

Maximum torque : 70 Ncm
Lighting mode : fiber optic
Lighting color : natural daylight
Cooling mode : external spray
Chuckling mode : push-button



VENTON M3 LED IMPLANTOLOGY MOTOR

High-end performance



Warranty :
2 years for the device
1 year for the contra-angle

Maximum torque	70 Ncm
Screen	Touch screen
Speed and torque display	Real time
Adjustable motor range (rpm)	100 - 40.000 (rpm)
Handpiece torque control	Yes
Coolant volume at 100%	> 90ml/min
Transmission ratio	1:1 20:1
Program settings	10 programs for individual adjustment
Patient files	Allows storage and export of data
PV	220V+/-10% - 50~60Hz



Motor

Large touch screen, new integrated pump for easy tube installation, automatic wire cutting mode
10 Customizable and free programming to meet different user needs

Foot control

Functions: Water Volume switch, program conversion, motor rotating direction control, stepless speed control




Micro-motor

Compact design, ideal balance and comfort
5-70 Ncm Torque range with 20:1 handpiece

Green ring contra-angle 20:1 LED

Maximum torque : 70 Ncm
Lighting mode : fiber optic
Lighting color : natural daylight
Cooling mode : external spray
Chuckling mode : push-button





Espace 14 it's...

ACCESS IMPLANTS

The Access implant system is manufactured in Germany and is integrated with the latest digital technologies (Exocad, 3Shape, Realguide, etc.)



Espace 14 it's also...



Biomaterials



Actilink
Plasma Activator



Instrumentations



Piezo Surgery
motors



G.B.R



Co2 / Diode lasers

TRAINING

Espace 14 has established a partnership with KCF Paris, a leading provider of dental education.



International Master's programs, tailored clinical coaching, short courses lasting 1 to 3 days on various topics (with or without official accreditation in France), and one-week clinical coaching programs in international clinics...

KCF offers training in oral implantology and dental innovation, tailored to the needs of professionals, with programs designed to help you deepen your skills and stay up to date with the latest techniques.

ESPACE 14

14 Place Jules Ferry
69006 Lyon - France

www.espace14-dental.com

Informations and contacts

Tel : 0033 (0)4 37 24 32 41
Tel : 0033 (0)6 28 68 07 02

contact@espace14-dental.com